

# Lovibond® Colour Measurement Tintometer® Group



## AF 710-3

Visual Grading for the Wesson Method



- Visual Grading of Oils & Fats
- Approved by the American Oil Chemists Society (AOCS)
- Measures According to AOCS Official Method Cc 13b-45

Part Number: 107103

### Visual Colorimeter for AOCS Official Method Cc 13b-45, the Wesson Method

The AF 710 is a visual colorimeter designed to optimise the use of Lovibond® glass filters in order to express colour in terms of the AOCS-Tintometer® Colour Scale. The AOCS-Tintometer® Colour Scale is a special red and yellow version of the Lovibond® RYBN scale using the AOCS-Tintometer® Red Scale (which is not the same as the Lovibond® standard red scale). The AF 710-3 accommodates the special flat-bottomed glass sample tubes specified in the AOCS Official Method Cc 13b-45, which are marked to indicate an oil column of ½", 1" and 5¼".

### International Standards Approval

The Lovibond® AF 710-3 is a visual instrument approved by the American Oil Chemists Society (AOCS) for colour grading of oils and fats.

### Easy Determination of Colour

You can easily measure the colour of the sample by comparing light transmitted through a column of oil with that transmitted through a series of red and yellow glasses calibrated in accordance with the AOCS-Tintometer® Colour Scale. The glasses are varied until a colour match is found for the light passing through the oil. You can then express this colour in units of red and yellow.

### No Excessive Ambient Light

The tube holder, which prevents light from entering the sides of the tubes, houses two tubes, one containing the sample and the other allowing an equal amount of light to reach the colour glasses.

## Industry

Chemical Industry | Food and Beverage Industry | Industries Others | Oil and Gas

## Application

Edible Oils and Fats | Household Products | Petroleum Oils and Waxes

## AF 710-3

The AF710-3 is a visual instrument approved by the American Oil Chemists Society (AOCS) for colour grading of oils and fats according to AOCS Official Method Cc 13b-45, the Wesson Method.

## Technical Data

<b>Applicable Standard</b>	AOCS Cc 8d-55 AOCS Cc 13b-45 AOCS Cc 13j-97
<b>Color Type</b>	Transparent
<b>Power Supply</b>	230 V / 50-60 Hz or 115 V / 50-60 Hz
<b>Portability</b>	Benchtop
<b>Touchscreen</b>	No

## Accessories

Title	Part Number
AF70-3 AOCS Racks	109080
Drip Tray	114199
Porcelain for AOCS	115600
Lamp 12v20w	122340

## Delivery Scope

AF 710-3 is supplied complete with a set of AOCS-Tintometer colour glasses, two glass sample tubes and a calibrated glass filter for regular conformance testing.

### Tintometer GmbH

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Germany

### The Tintometer Limited

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
UK

### Tintometer China

Room 1001, China Life Tower  
16 Chaoyangmenwai Avenue,  
Beijing, 100020  
Tel.: +86 10 85251111 App. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
China

### Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,  
Lebuh Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Malaysia

### Tintometer Brazil

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@tintometer.com.br  
www.lovibond.com.br  
Brazil

### Tintometer Inc.

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.com  
USA

### Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
India

### Tintometer Spain

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
Spain

Technical changes without notice  
Printed in Germany

Lovibond® and Tintometer® are Trademarks of the Tintometer Group of Companies