

315 The new biological microscope can observe cells 7-inch large-screen digital microscope 2000x high magnification microscope










No. 1 lens




No. 2 lens





- | | | |
|----------------------------------|--|-----------------------------------|
| 1 7 inch ips display | 7 Focus adjustment roller | 13 Display fixing nut |
| 2 Charging indicator | 8 Bracket adjustment knob | 14 18650 battery slot |
| 3 Micro TF card interface | 9 Object Stage | 15 Image adjustment wheel |
| 4 Micro USB socket | 10 Bracket LED adjustment wheel | 16 LED adjusting wheel |
| 5 USB power output | 11 Bracket light source | 17 Bracket power interface |
| 6 Reset hole | 12 Removable lens | |



Real 2000X effect

False 1000X effect

others

ours





1



8



2



3



4



5



7



6

1. Microscope x1
2. Low multiple lens x1
3. High multiple lens x1
4. Manual x1
5. Slide x2
6. Data cable x2
7. Cleaning cloth x1
8. Packing box x1